



Product Code	34408
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No

RDM (NSIS)	2334266
CND	R01020101
EAN/UPC	5707480023172
PARAF	988194650
GMDN	45036

## Description

### AMBU® AURAONCE™ DISPOSABLE LARYNGEAL MASKS

- Laryngeal Mask N 3
- Colour: green
- Patient weight: 30-50 kg
- Maximum cuff inflation volume: 20 ml

AuraOnce has a special curve that carefully replicates natural human anatomy.

Molded directly into the tube, the unique shape makes correct insertion easy without abrading the upper airway.

AuraOnce's curve also ensures that the patient's head remains in a natural supine position when the mask is in use without extra stress on the upper jaw. Moreover, Auraonce is specially designed to give the airway tube the flexibility needed to adapt to individual anatomical variances and wide range of head positions.

AuraOnce has a D-shaped airway tube to give a firm and ergonomical grip during insertion. In short, positioning is fast and accurate, with least possible internal pressure.

The airway tube is flexible at the cuff and rigid at the connector for easy, atraumatic insertion and removal.

Cuff and airway tube moulded as single unit with built-in, anatomically correct curve.

Smooth sides without ridges or fins that can scratch delicate tissue.

Reinforced tip will not bend during insertion so positioning is always correct.

Ultra thin pilot balloon with universal check valve provides precise tactile indication of degree of inflation.

Reinforced tip will not bend during insertion so positioning is always correct.

Universal 15 mm connector (ISO).

## Technical Specifications

Maximum intra-cuff pressure: 60 cm H<sub>2</sub>O

Airway connector: 15 mm ISO 5356-1

Tube/cuff material: PVC (latex and phthalates free)

Complete range for neonatal (AuraGain only), pediatric and adult patients

Packaged sterile and ready for use

Colour coded pouch with directions for use