



| | |
|---------------|----------------|
| Product Code | 34418 |
| Unit of sale | 1 pc |
| Minimum order | 1 |
| Type | Medical device |
| Class | I |

| | |
|--------|-----|
| UK-REP | Yes |
|--------|-----|

| | |
|------------|---------------|
| RDM (NSIS) | 2789843 |
| CND | Z12021003 |
| EAN/UPC | 4053755127481 |
| UMDNS | 15-076 |

Description

HEINE FIBRE OPTIC BLADES

- HEINE Code: F-000.22.101
- Blade type: MAC 1
- Overall length: 93 mm
- Distal width: 11 mm

The "Classic+" range of F.O. blades with HEINE HiLite®

Fibre bundle with over 6,500 individual micro-fibres offer improved light transmission and longer life. Together with the innovative HEINE LEDHQ light in the handle or the traditional XHL Xenon Halogen Bulb, the exclusive HEINE HiLite® Fibre Optic bundle with up to 4.3 mm of Ø offers 3X or 40% brighter light as with other conventional halogen bulbs.

5-years guarantee or 4,000 standard autoclave cycles of 134°C/5 min. with more than 1,000 lux on all blades with HEINE HiLite® fibre bundle.

Long lasting bulbs

HEINE XHL® Xenon Halogen bulbs incorporate the latest Xenon Halogen Technology to give you up to 40% more light! XHL 3.5 V bulbs last up to twice as long as 2.5 V bulbs and give up to three times the light of conventional 2.5 V Halogen laryngoscope bulbs. HEINE LEDHQ handles offer even three times the light intensity of the 3.5V Xenon Halogen bulbs for an outstanding and constant bright illumination.

Long lasting 3.5 V rechargeable batteries

HEINE 3.5 V systems (LEDHQ or XHL) use a long-lasting 3.5 V lithium rechargeable battery. HEINE F.O. laryngoscope blades and handles are compatible with the ISO 7376 (green standard).