



Product Code	44300
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
--------	-----

RDM (NSIS)	2225488
CND	V08020101
EAN/UPC	8023279443004

Description

MAYA GYNAECOLOGICAL CHAIRS - blue

Made of tubular furniture with remarkable performance characteristics that guarantee stability and sturdiness. The very compact base of the structure allows the operator to get closer to the patient, making it very practical and comfortable to use. The structure has excellent maneuverability determined by the use of 3 electric actuators controlled by a push-button. Central column allows to adjust the structure in height, from 50 to 100 cm. Adjustable backrest by electric actuator to reach +75° and -10°; the leg cushion reaches -90° and +5°. The upper surface and the leg loops are in high density foam covered in coated fabric (skai, imitation leather) fireproof class 11M, washable with common detergents on the market. Rubber density guarantees comfort and proper support to the patient during the examination. Supplied with a paper roll holder and 4 adjustable feet that guarantee greater stability.

Size: 182x68xh 57-98 cm.

Weight: 80 kg.

Max. load: 250 kg.

Made in Italy.