



Product Code	44320
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	2843600
CND	V08020199
EAN/UPC	8023279443202

Description

ELECTRIC DIALYSIS CHAIR

Comfortable chair ideal for patients who have to undergo long dialysis sessions. Steel frame, plastic seat and flame-retardant PVC leather. Mattress thickness: 6 cm.

Leg rest, backrest and height can be electronically adjusted by hand controller. Armrests are manually adjustable: right/left, front/back, for patient's comfort. Integrated adjustable ABS table board, with antiskid convex edge and holes to place cups and cell phone.

The chair is also equipped with CPR quick release system and 4 silent castors Ø 76 mm with brake.

Multilingual instructions:

Max load: 240 kg.

Size: 200x58x h 55-70 cm.

Weight: 100 kg.

Voltage: 100-240 V - 50/60 Hz.

Backrest size: 77x58 cm

Backrest lifting: 8° to 70°

Seat size: 47x58 cm

Footrest size: 47x23 cm

Footrest can extend chair length by 25 cm

Armrest size: 50x18 cm

Armrest extension: 0° to 270°

Leg rest size: 40x54 cm

Leg rest lifting: -60° to 0°