

EXAMINATION COUCH load 120 kg - white painted, black mattress



Product Code	44452
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2537737
CND	V080603
EAN/UPC	8023279444520
PARAF	989629290
GTIN	08023279444520
GMDN	38458

Description

VALUE EXAMINATION COUCH - white painted - black

Suitable for medical practices or small surgeries.

Top is made with wooden panel, foam rubber and covered with washable, fireproof simulated leather Class 1.

Two sections: one fixed 131x62xh 6 cm and one adjustable (backrest) 51x62xh 6 cm.

Couch size: 182x62xh 79 cm.

Weight 27 kg.

Maximum load: 120 kg.

Delivered in kit form: 173x74xh 16 cm.

Easy to assemble.

Multilingual manual (GB, FR, IT, ES, PT, DE).

Made in Italy.