

CHROMED PLATED STEEL EXAMINATION COUCH



Product Code	44461
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No
RDM (NSIS)	521354
CND	Y181207
EAN/UPC	8023279444612
PARAF	988194623

Description

Chrome plated examination couch

Structure made of chrome plated steel tube Ø 25x25 mm.

Backrest adjustable by rack.

Padding and covering made of washable, flame retardant, plastic material. Rubber feet.

Supplied with a paper roll holder, which can be placed either in the front or in the back part of the examination couch.

Fireproof padding and covering, Class 1 IM.

Dimension: 177 x 60 x h 75 cm.

Weight: 25 kg

Load 150 Kg

Delivered disassembled.

Made in Italy.