

STAINLESS STEEL EXAMINATION COUCH



| | |
|---------------|----------------|
| Product Code | 44463 |
| Unit of sale | 1 pc |
| Minimum order | 1 |
| Type | Medical device |
| Class | I |

| | |
|--------|----|
| UK-REP | No |
| CH-REP | No |

| | |
|------------|---------------|
| RDM (NSIS) | 470792 |
| CND | Y181207 |
| EAN/UPC | 8023279444636 |

Description

Stainless steel examination couch

Structure made of stainless steel tube Ø 25x25 mm.

Backrest adjustable by rack.

Padding and covering made of washable, flame retardant, plastic material. Rubber feet.

Supplied with a paper roll holder, which can be placed either in the front or in the back part of the examination couch.

Fireproof padding and covering, Class 1 IM.

Dimension: 177 x 60 x h 75 cm.

Weight: 25 kg

Load 150 Kg

Delivered disassembled.

Made in Italy.