

# STAINLESS STEEL EXAMINATION COUCH WITH 3 DRAWERS



Product Code	44467
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	2635498
CND	V0899
EAN/UPC	8023279444674
GMDN	38459

## Description

### STAINLESS STEEL EXAMINATION COUCH WITH 3 DRAWERS

Structure made of enamelled stainless steel tube Ø 30x30 mm.

Adjustable backrest by rack.

Padding and covering made of washable, flame retardant, plastic material.

Both equipped with three drawers to be placed under the fixed part.

Size: 177x60xh 80 cm.

Drawer internal size: 38x45x h 10 cm.

Single drawer load capacity: 10 kg.

Packing size 1 box 74x21x202 cm + 1 box 77x51x55 cm.

Class 1 IM.

Multilingual manual and instructions: GB, FR, IT, ES PT, DE, NL, DK, SE, PL, CZ, HU, SI, HR, RO, BG, GR, Arabic.

Delivered disassembled.

Made in Italy.