

# ELITE EXAMINATION COUCH - blue



Product Code	44500
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	1975772
CND	V0899
EAN/UPC	8023279445008

## Description

### ELITE EXAMINATION COUCH - blue

Stylish, comfortable and innovative examination couch in blue color with lateral white band.

Realized with high quality furniture tubing that guarantees stability and robustness.  
The couch is made with rectangular and rounded tubing, Uni 7947/81 and 8913, DIN 2394.  
Fireproof rubber padding covered with class 1 IM fireproof woven imitation leather.  
Adjustable head section through 2 cogwheels.  
The roll holder, included, can be placed in the feet or in the head side.  
4 adjustable feet.

Size: 193x68xh 78 cm.  
Padding Size: 190x68xh 7 cm.  
Max. load 200 kg.  
Weight 70 kg.  
Disassembled, very simple to assemble.

Made in Italy.