


LORD HEIGHT ADJUSTABLE EXAMINATION COUCH with TR/RTR - water green



Product Code	44510
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	1793166
CND	V0899
EAN/UPC	8023279445107

Description

LORD VARIABLE HEIGHT EXAMINATION COUCH - water green

Stable and sturdy examination couch made of high quality tubular for furniture. Very practical and user friendly. Height adjustable thanks to the two electrical columns that can move simultaneously or separately (allowing the Trendelenburg/Reverse Trendelenburg position) bby a push-button panel. Under the top frame, an ergonomic lever activates a gas spring that reclines the cushions. High density, rubber padded, fire retardant and Class 1 IM top, covered in coated fabric (vinyl, synthetic leather, flame retardant Class 1 IM). It can be cleaned with common detergents. Complies with: EU 93/42 CEE.

Maximum load: 200 kg.
Supplied with mouth/nose hole plug, paper roll holder and 4 adjustable feet.

Overall size:
194x68xh 52/94 cm
Padded top length: 190 cm
Weight: 80 kg

Made in Italy.