

ELECTRIC HEIGHT ADJUSTABLE TREATMENT TABLE with footbar - blue



| | |
|---------------|----------------|
| Product Code | 44520 |
| Unit of sale | 1 pc |
| Minimum order | 1 |
| Type | Medical device |
| Class | I |

| | |
|--------|-----|
| UK-REP | No |
| CH-REP | Yes |

| | |
|------------|----------------|
| RDM (NSIS) | 2537740 |
| CND | V08060199 |
| EAN/UPC | 8023279445206 |
| GTIN | 08023279445206 |
| GMDN | 38447 |

Description

Electric table - footbar - blue

Sturdy and highly stable examination and treatment table.
Adjustable height by electric footbar.
Adjustable head-section by gas pump with mouth-nose hole.
Sturdy tubular steel frame painted with non-toxic epoxy powder, RAL 9010.
Chromed footbar.

Intended for professional use in hospitals, physiotherapy, treatments.

Supplied with 6 cm thickness mattress and mouth-nose hole plug 23x11 cm included.

Fireproof padding.
Available in 12 colours.
Class I, Medical Devices, complies with MDR 2017/745.
Maximum overall dimensions 195 x 59,5 x h 54-71 cm.
Weight: 65 kg.
Lifting capacity: 135 kg.
Multilingual manual: GB, FR, IT, ES, PT, DE, GR.

Made in Italy.