

# ELECTRIC HEIGHT ADJUSTABLE TREATMENT TABLE with footbar - black



Product Code	44522
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	Yes

RDM (NSIS)	2537740
CND	V08060199
EAN/UPC	8023279445220
GTIN	08023279445220
GMDN	38447

## Description

Electric table - footbar - black

Sturdy and highly stable examination and treatment table.  
Adjustable height by electric footbar.  
Adjustable head-section by gas pump with mouth-nose hole.  
Sturdy tubular steel frame painted with non-toxic epoxy powder, RAL 9010.  
Chromed footbar.

Intended for professional use in hospitals, physiotherapy, treatments.

Supplied with 6 cm thickness mattress and mouth-nose hole plug 23x11 cm included.

Fireproof padding.  
Available in 12 colours.  
Class I, Medical Devices, complies with MDR 2017/745.  
Maximum overall dimensions 195 x 59,5 x h 54-71 cm.  
Weight: 65 kg.  
Lifting capacity: 135 kg.  
Multilingual manual: GB, FR, IT, ES, PT, DE, GR.

Made in Italy.