

4-SECTIONS HEIGHT ADJUST. TREATMENT COUCH - blue



Product Code	44531
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
RDM (NSIS)	2346026
CND	V0899
EAN/UPC	8023279445312

Description

4-SECTION VARIABLE HEIGHT TREATMENT COUCH - blue

Couch made of high quality tubular that makes it very sturdy.

The cushions can be reclined by a gas spring activated by an ergonomic lever under the top frame.

Head section movable to +70° and to -35°.

The two electrical columns allow to adjust height and reach TR/RTR position very quickly by pedals, push-button or two perimeter bars.

High density, rubber padded, fire retardant and Class 1 IM top, covered in flame retardant Class 1 IM coated fabric (synthetic leather).

Washable and disinfectable with common detergents.

Overall size: 195x68xh 48/89 cm. Weight: 80 kg.

Complies with: EU 93/42 CEE.

Maximum load: 200 kg.

Supplied with mouth/nose hole plug, paper roll holder, 2 adjustable armrests and 4 adjustable feet.

Made in Italy.