

ROYAL ADJUSTABLE EXAMINATION COUCH - stainless steel



Product Code	44541
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	2879102
CND	V080699
EAN/UPC	8023279445411

Description

ROYAL ADJUSTABLE EXAMINATION COUCH - blue

Couch frame is made of painted powder coated steel with strong corrosion resistance: this structure allows to maintain more stability, making patients feel safe.

Couch surface is made of leather and hand-made cotton mattress, which make it soft and breathable.

Easy to clean and disinfect.

It is equipped with manual air levers with two gas springs to adjust the back and leg plate angle, one hand controller to adjust height.

The brake system is composed by four brake pedals connected to four castors Ø 76 mm.

360° spinning castors made of strong, durable and non-slip polyurethane PU.

Size: 196x68xh 47-91 cm. Weight: 68 kg.

Safe working load: 160 kg.

Voltage: AC 220 V, 50 Hz.