



Product Code	44690
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No
RDM (NSIS)	35032
CND	Z110702
EAN/UPC	8023279446906

Description

GIMA X-ray illuminators have been designed according to clinical recommendations for an extra white light to give an accurate and precise reading of any detail.

Made of enamelled sheet steel.

Unbreakable white monochromatic diffuser, stainless steel gripping system.

Made in Italy.

- Homogeneous fluorescent cold light
- Brightness = 1,700 cd/m²
- Power supply = 230 V - 50/60 Hz

- External dimensions

- Length: 43 cm
- Height: 46 cm
- Depth: 14 cm
- Weight: 6 Kg
- Type: 1 panel

Made in Italy.