



| | |
|---------------|----------------|
| Product Code | 44692 |
| Unit of sale | 1 pc |
| Minimum order | 1 |
| Type | Medical device |
| Class | I |

| | |
|------------|---------------|
| UK-REP | No |
| CH-REP | No |
| RDM (NSIS) | 35032 |
| CND | Z110702 |
| EAN/UPC | 8023279446920 |

Description

GIMA X-ray illuminators have been designed according to clinical recommendations for an extra white light to give an accurate and precise reading of any detail.

Made of enamelled sheet steel.

Unbreakable white monochromatic diffuser, stainless steel gripping system.

On/off switch for separate use of each panel.

Made in Italy.

- Homogeneous fluorescent cold light
- Brightness = 1.700 cd/m²
- Power supply = 230 V - 50/60 Hz

- External dimensions
 - Length: 83 cm
 - Height: 46 cm
 - Depth: 14 cm
- Weight: 10 Kg
- Type: 2 panels

Made in Italy.