

MASTER OVERBED TABLE



Product Code	27447
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

EAN/UPC	8023279274479
---------	---------------

Description

MASTER OVERBED TABLE

High quality overbed table with adjustable height by gas spring system.
Easy to clean table top made of ABS material, equipped with cup holder and high edges to protect objects from falling down.
Anti-collision and durable round corner design.
Bottom frame is sprayed with carbon steel.
4 revolving rubber castors Ø 50 mm, 2 with brake.

Size: 890x440xh 765-1050 mm.

Height adjustment (765-1050 mm)
Extra large abs table top with cup socket
Air spring for height adjustment