



Product Code	44701
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	2588211
CND	Z110702
EAN/UPC	8023279447019

Description

WALL LED SLIM VIEWER 42 x 72 cm

- External size
- Length: 82 cm
- Height: 50,5 cm
- Depth: 3.8 cm
- Weight: 8 Kg
- Switches: 1
- Types: 2 pannels

LED TECHNOLOGY SLIM X-RAY ILLUMINATORS

The new generation of LED X-ray illuminators with colour temperature of 7,200°K for superior viewing and reading of radiography. Complete elimination of light flickering and 5,000 cd/mq luminance for more uniform light diffusion. Frame in antiscratch anodized aluminium. Film gripping system incorporated.

Dimmer for luminance adjustment from 30% to 100%.

- FAS: automatic switching on/off at the positioning/removal of the film.

Economic advantages:

- LED life 80,000 hours, with low maintenance costs
- low energy consumption compliant with European law

Made in Italy.

Technical Specifications

Colour temperature: Cool white 7,200°K
 Brightness: 5,000 cd/m²
 Power supply: 110-230 V - 50/60 Hz
 LED life: 80,000 hours
 Conformity to Low Voltage Directive 2014/35/CE
 Conformity to Electro Magnetic compatibility 2014/30/CE
 Conformity to EEC directive 2007/47 medical Device Class I
 All components UL, DVE, IMQ approved