



|               |                |
|---------------|----------------|
| Product Code  | 44705          |
| Unit of sale  | 1 pc           |
| Minimum order | 1              |
| Type          | Medical device |
| Class         | I              |

|            |               |
|------------|---------------|
| UK-REP     | No            |
| CH-REP     | No            |
| RDM (NSIS) | 1040806       |
| CND        | Z110702       |
| EAN/UPC    | 8023279447057 |

## Description

DESKTOP LED SLIM VIEWER 42 x 36 cm x1

- External size
  - Length: 46 cm x1
  - Height: 50,5 cm
  - Depth: 3.8 cm
- Weight: 6 Kg
- Switches: 1
- Types: 1 pannel

### LED TECHNOLOGY SLIM X-RAY ILLUMINATORS

The new generation of LED X-ray illuminators with colour temperature of 7,200°K for superior viewing and reading of radiography. Complete elimination of light flickering and 5,000 cd/mq luminance for more uniform light diffusion. Frame in antiscratch anodized aluminium. Film gripping system incorporated.

Dimmer for luminance adjustment from 30% to 100%.

- FAS: automatic switching on/off at the positioning/removal of the film.

Economic advantages:

- LED life 80,000 hours, with low maintenance costs
- low energy consumption compliant with European law

Fully adjustable desktop model:

Position of panel can be totally adjusted:

- height adjustment 13 cm
- inclination 35° (+25°/-5°)
- vertical or horizontal positioning
- horizontal rotation 270°

Made in Italy.

## Technical Specifications

Colour temperature: Cool white 7,200°K  
 Brightness: 5,000 cd/m<sup>2</sup>  
 Power supply: 110-230 V - 50/60 Hz  
 LED life: 80,000 hours  
 Conformity to Low Voltage Directive 2014/35/CE  
 Conformity to Electro Magnetic compatibility 2014/30/CE  
 Conformity to EEC directive 2007/47 medical Device Class I  
 All components UL, DVE, IMQ approved