

DESKTOP ULTRA SLIM LED LIGHT BOX 42x36 cm - double



Product Code	44706
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I
UK-REP	No
CH-REP	No
RDM (NSIS)	1040826
CND	Z110702
EAN/UPC	8023279447064

Description

DESKTOP LED SLIM VIEWER 42 x 36 cm x2

- External size
 - Length: 46 cm x2
 - Height: 50.5 cm
 - Depth: 3.8 cm
- Weight: 20 Kg
- Switches: 2
- Types: 2 pannels

LED TECHNOLOGY SLIM X-RAY ILLUMINATORS

The new generation of LED X-ray illuminators with colour temperature of 7,200°K for superior viewing and reading of radiography. Complete elimination of light flickering and 5,000 cd/mq luminance for more uniform light diffusion. Frame in antiscratch anodized aluminium. Film gripping system incorporated.

Dimmer for luminance adjustment from 30% to 100%.
- FAS: automatic switching on/off at the positioning/removal of the film.

Economic advantages:
- LED life 80,000 hours, with low maintenance costs
- low energy consumption compliant with European law

Fully adjustable desktop model:
Position of panel can be totally adjusted:
- height adjustment 13 cm
- inclination 35° (+25°/-5°)
- horizontal rotation 10°

Made in Italy.

Technical Specifications

Colour temperature: Cool white 7,200°K
Brightness: 5,000 cd/m²
Power supply: 110-230 V - 50/60 Hz
LED life: 80,000 hours
Conformity to Low Voltage Directive 2014/35/CE
Conformity to Electro Magnetic compatibility 2014/30/CE
Conformity to EEC directive 2007/47 medical Device Class I
All components UL, DVE, IMQ approved