

# DESKTOP ULTRA SLIM LED LIGHT BOX 42x36 cm - triple



Product Code	44707
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	1040846
CND	Z110702
EAN/UPC	8023279447071

## Description

DESKTOP LED SLIM VIEWER 42 x 36 cm x3

- External size
  - Length: 46 cm x3
  - Height: 50,5 cm
  - Depth: 3.8 cm
- Weight: 26 Kg
- Switches: 3
- Types: 3 pannels

### LED TECHNOLOGY SLIM X-RAY ILLUMINATORS

The new generation of LED X-ray illuminators with colour temperature of 7,200°K for superior viewing and reading of radiography. Complete elimination of light flickering and 5,000 cd/mq luminance for more uniform light diffusion. Frame in antiscratch anodized aluminium. Film gripping system incorporated.

Dimmer for luminance adjustment from 30% to 100%.  
- FAS: automatic switching on/off at the positioning/removal of the film.

Economic advantages:  
- LED life 80,000 hours, with low maintenance costs  
- low energy consumption compliant with European law

Fully adjustable desktop model:  
Position of panel can be totally adjusted:  
- height adjustment 13 cm  
- inclination 35° (+25°/-5°)  
- horizontal rotation 10°

Made in Italy.

## Technical Specifications

Colour temperature: Cool white 7,200°K  
Brightness: 5,000 cd/m<sup>2</sup>  
Power supply: 110-230 V - 50/60 Hz  
LED life: 80,000 hours  
Conformity to Low Voltage Directive 2014/35/CE  
Conformity to Electro Magnetic compatibility 2014/30/CE  
Conformity to EEC directive 2007/47 medical Device Class I  
All components UL, DVE, IMQ approved