



Product Code	44711
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No
RDM (NSIS)	2149603
CND	Z12120199
EAN/UPC	8023279447118

## Description

New generation of ultra slim LED optometric charts.  
Made in Italy.

- Monoyer - 3 m
- Viewing area: 25.5 x 63 cm
- External size: 34 x 72 x 3.8 cm

## Technical Specifications

Colour temperature: Cool white 7,200°K  
 Brightness: 5,000 cd/m<sup>2</sup>  
 Power supply: 110-230 V - 50/60 Hz  
 LED life: 80,000 hours  
 Conformity to Low Voltage Directive 2014/35/CE  
 Conformity to Electro Magnetic compatibility 2014/30/CE  
 Conformity to EEC directive 2007/47 medical Device Class I  
 All components UL, DVE, IMQ approved