



Product Code	44712
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No
RDM (NSIS)	1346229
CND	Z12120199
EAN/UPC	8023279447125

Description

New generation of ultra slim LED optometric charts.
Made in Italy.

- Mixed decimal - 3m
- Viewing area: 25.5 x 63 cm
- External size: 34 x 72 x 3.8 cm

Technical Specifications

Colour temperature: Cool white 7,200°K
 Brightness: 5,000 cd/m²
 Power supply: 110-230 V - 50/60 Hz
 LED life: 80,000 hours
 Conformity to Low Voltage Directive 2014/35/CE
 Conformity to Electro Magnetic compatibility 2014/30/CE
 Conformity to EEC directive 2007/47 medical Device Class I
 All components UL, DVE, IMQ approved