

# 3 JOINTS VARIABLE HEIGHT PATIENT BED - 2 cranks - castors



Product Code	44734
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	1352082
CND	Y181299
EAN/UPC	8023279447347

## Description

3 JOINTS BED - 2 cranks - variable height - castors

- 3 joints/4 sections beds adjustable by means of 2 independent cranks
- height adjustable by hydraulic foot pump
- all welded main frame constructed from rectangular steel frame
- bed ends made in epoxy finish steel tube
- legs fitted with 4 revolving rubber/steel castors Ø 25mm with independent brake

## Technical Specifications

Size of platform: 900 x 2,000 mm  
Platform height mm: 400/800  
Head/foot board height mm: 460  
External measurement: 910 x 2,060 mm  
Movement: 2 cranks  
Max patient load kg: 135