

# 3 JOINTS PATIENT BED - 2 cranks - castors



|               |                |
|---------------|----------------|
| Product Code  | 44740          |
| Unit of sale  | 1 pc           |
| Minimum order | 1              |
| Type          | Medical device |
| Class         | I              |

|        |    |
|--------|----|
| UK-REP | No |
| CH-REP | No |

|            |               |
|------------|---------------|
| RDM (NSIS) | 2321126       |
| CND        | V080699       |
| EAN/UPC    | 8023279447408 |

## Description

### 3 JOINTS PATIENT BED - 2 cranks - castors

Manual patient beds with 3 joints adjustable by means of 2 cranks.  
PP head and foot board equipped with bumpers wheel Ø 65 mm, preventing collision.  
PP bed platform, divided in 4 sections with integrated mattress retainer.  
Plastic parts and epoxy painting.

#### Guardrail

Aluminium alloy side rail, 6 staff, collapsed type with safe lock

#### Castors

4 Ø 12.7 cm castors

#### Mattress retainer

Platform integrated mattress retainer on both sides (back and leg section), with a boarder of 33 mm high

#### IV Pole

Stainless steel IV pole, bearing 15 kg

## Technical Specifications

Size: 2130x995 mm  
Height: 500 mm  
Mattress platform: 1925x900 mm  
Max. patient load: 200 kg  
Movement: 2 cranks  
Back-rest adjustment: 0-70° (± 5°)  
Knee-rest adjustment: 0-38° (± 5°)

## Standard accessories

2 foldable side rails  
IV pole 4 hooks  
4 IV pole provisions  
2 drainage hooks