

3 JOINTS PATIENT BED - 2 cranks - castors



Product Code	44740
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	2475745
CND	V080602
EAN/UPC	8023279447408

Description

3 JOINTS PATIENT BED - 2 cranks - castors

Manual patient beds with 3 joints adjustable by means of 2 cranks.
PP head and foot board equipped with bumpers wheel Ø 65 mm, preventing collision.
PP bed platform, divided in 4 sections with integrated mattress retainer.
Plastic parts and epoxy painting.

Guardrail

Aluminium alloy side rail, 6 staff, collapsed type with safe lock

Castors

4 Ø 12.7 cm castors

Mattress retainer

Platform integrated mattress retainer on both sides (back and leg section), with a boarder of 33 mm high

IV Pole

Stainless steel IV pole, bearing 15 kg

Technical Specifications

Size: 2130x995 mm
Height: 500 mm
Mattress platform: 1925x900 mm
Max. patient load: 200 kg
Movement: 2 cranks
Back-rest adjustment: 0-70° (± 5°)
Knee-rest adjustment: 0-38° (± 5°)

Standard accessories

2 foldable side rails
IV pole 4 hooks
4 IV pole provisions
2 drainage hooks