

4 JOINTS VARIABLE HEIGHT PATIENT BED - 3 cranks - castors



| | |
|---------------|----------------|
| Product Code | 44742 |
| Unit of sale | 1 pc |
| Minimum order | 1 |
| Type | Medical device |
| Class | I |

| | |
|------------|---------------|
| UK-REP | No |
| CH-REP | No |
| RDM (NSIS) | 2321126 |
| CND | V080699 |
| EAN/UPC | 8023279447422 |

Description

4 JOINTS PATIENT BED - 3 cranks - variable height - castors

Manual patient beds with 4 joints adjustable by means of 3 cranks.
PP head and foot board equipped with bumpers wheel Ø 65 mm, preventing collision.
PP bed platform, divided in 4 sections with integrated mattress retainer.
Plastic parts and epoxy painting.

Guardrail
Aluminium alloy side rail, 6 staff, collapsed type with safe lock

Castors
4 Ø 12.7 cm castors

Mattress retainer
Platform integrated mattress retainer on both sides (back and leg section), with a boarder of 33 mm high

IV Pole
Stainless steel IV pole, bearing 15 kg

Technical Specifications

Size: 2130x995 mm
Height: adjustable 420-730 mm
Mattress platform: 1925x900 mm
Max. patient load: 200 kg
Movement: 3 cranks
Back-rest adjustment: 0-70° (± 5°)
Knee-rest adjustment: 0-38° (± 5°)

Standard accessories

2 foldable side rails
IV pole 4 hooks
4 IV pole provisions
2 drainage hooks