



Product Code	44744
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No
RDM (NSIS)	2549027
CND	V080699
EAN/UPC	8023279447446

Description

ELECTRIC PATIENT BED - adjustable height

Electric patient beds with integrated PP bed platform and electric motor to facilitate different movements of bed sections. Hand controller allows nurses to control standard functions. Both models have an American motor with UL/ROHS/EN standard. Plastic parts and epoxy painting. PP head and foot board equipped with four bumpers wheels Ø 65 mm.

Mattress retainer

Platform integrated mattress retainer on both sides (back and leg section) with a boarder of 33 mm high

Guardrail

Aluminium alloy foldable side rail, 6 staff, collapsed type with safe lock

IV Pole

Stainless steel IV stand, bearing 15 kg

Castors

4 Ø 12.7 cm (5") anti-winding hard shell covered castors

Hand controller

It allows to adjust bed sections

Bumper wheel

Equipped with 4 bumpers Ø 65 mm to protect bed and wall

Technical Specifications

Size: 208.5x99.5 cm
 Height adjustment: from 42 to 73 cm
 Mattress platform: 192.5x90 cm
 Safe working load: 200 kg
 Back-rest adjustment: 0-75° (± 5°)
 Knee-rest adjustment: 0-46° (± 5°)
 Trendelenburg: -
 Reverse Trendelenburg: -
 Angle indicator: -
 Auto-regression: -
 Position line: -
 Motor life: 20.000 movements
 Motor noise: under 50 dB
 Voltage: 100-240 V - 50/60 Hz

Standard accessories

- 1 hand controller
- 1 set foldable side rails
- 1 IV stand with 4 hooks
- 4 IV pole provisions
- 2 drainage hooks