



Product Code	44770
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	2315777
CND	Y093312
EAN/UPC	8023279447705

Description

SHOWER TROLLEY - hydraulic

Shower stretcher can be used for both transport and showering. Ideal solution for treatment of patients in assisted baths of hospital or retirement home. Rubber castors Ø 12.5 cm with individual brake and one with directional lock to facilitate transport in long linear path. Automatic progressive tilting to allow a gradual water drainage. Safety side rails with one-hand folding/locking system. Supplied with PCM foam padded waterproof mattress and pillow made of fire-resistant auto-extinguishing Class 2 PVC, extendable and flexible drain hose (Ø 40 mm - length 1.5 m).

Technical Specifications

Max working load: 190 kg
Weight: 93 kg
Size: 2,000 x 745 x h 740//1080 mm
Mattress height: 20 cm
Tilt for water discharge: 7°
Discharge hose length: 150 cm
Collapsible head end rail: yes
Directional system: yes