



Product Code	44775
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	2549042
CND	Z12011202
EAN/UPC	8023279447750

## Description

### ELECTRIC OPERATING TABLE with TR and RTR

Electric operating table for general surgical procedures, mainly made of 304 stainless steel. It is driven by low-voltage DC motor with safe and stable performance. Table movements can be controlled by means of an easy-to-use hand controller. Super wide table top, super long sliding distance. Leg rest with integrated gas spring for easy control. Suitable for X-ray and C-arm surgery thanks to an extra-long horizontal sliding. Max. load: 200 kg

**Head and leg rest**  
Detachable double head and leg rest for surgical needs

**Hand controller**  
All functions can be controlled by hand controller and can be locked

**Kidney bridge**  
Built-in kidney bridge

**Five section table top**  
Head, back, shoulder hip and leg section, each section with X-ray cassette layer to allow X-ray use

## Technical Specifications

Size: 210x50xh 80-100 cm  
 Max. load: 200 kg  
 Trendelenburg: 0-25°  
 Reverse Trendelenburg: 0-22°  
 Lateral tilt (left/right): 0-20°  
 Head rest (up/down): 0-45°/0-90°  
 Back plate (up/down): 0-75°/0-20°  
 Foot plate (up/down): 0-15°/0-90°  
 Foot plate horizontal positioning: 0-90°  
 Horizontal sliding: 30 cm  
 Kidney bridge elevation: 11 cm  
 Voltage: 110-220 V, 50/60 Hz  
 Power capacity: 1.0 Kw

## Standard accessories

Matress

- 1 anesthesia screen frame
- 2 shoulder holders with cushion
- 2 arm holders with cushion
- 2 leg supports with cushion
- 2 foot supports with cushion
- 1 kidney bridge handle
- 1 IV stand