

ADJUSTABLE HEIGHT PATIENT TROLLEY



Product Code	44780
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No
RDM (NSIS)	2475887
CND	V08050101
EAN/UPC	8023279447804

Description

HEIGHT ADJUSTABLE PATIENT TROLLEY

Height adjustable patient trolley made of metal.
Both bed frame and chassis have 30x50 mm square tube to make stretcher stronger and more stable.
Bed platform and guardrail are made of high quality PP material.
Back lift makes patient more comfortable.
Equipped with strips to prevent patients from possibly sliding off the stretcher and 2 cm thickness mattress.
Directional wheel and central locking pedals.

Multilingual manual: GB, FR, IT, ES, PT, DE, GR, RO.

Directional central wheel to facilitate direction control.

Four double side castors Ø 15 cm lockable by pedal.

Technical Specifications

Size: 1900x640 mm
Height adjustment: 550-850 mm
Back-plate adjustment: 0-70° (± 10°)
Safety load: 170 Kg.

Standard accessories

1 set IV pole
1 transfer mattress
2 IV pole provision
2 strips
1 oxygen bottle holder for 2 l bottle
1 set ABS folding side rails