

HYDRAULIC ADJUSTABLE HEIGHT PATIENT TROLLEY with TR and RTR



Product Code	44782
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	2475889
CND	V08050101
EAN/UPC	8023279447828

Description

HYDRAULIC HEIGHT ADJUSTABLE PATIENT TROLLEY with TR and RTR

Metal patient trolley equipped with foldable aluminium alloy guardrails and two hydraulic pumps which allow height adjustment, Trendelenburg and reverse Trendelenburg.
Four double castors Ø 20 cm, lockable by pedals on both sides. Waterproof PU mattress.
Four detachable bumper wheels Ø 74 mm to protect bed and patients during transportation.
Stainless steel IV pole with plastic hooks.
X-ray translucent bed platform.
Central wheel to facilitate direction control.

Multilingual manual: GB, FR, IT, ES, PT, DE, GR, RO.

Technical Specifications

Size: 2100x850 mm
Height adjustment: 630-930 mm
Back-plate adjustment: 0-75° ($\pm 10^\circ$)
Safety load: 220 Kg.
Knee-rest adjustment: 0-45°
Trendelenburg: 0-15°
Reverse Trendelenburg: 0-15°

Standard accessories

1 set emergency push handle
4 IV pole provisions
1 set folding side rails
1 mattress
1 set IV pole