


BARIATRIC MATTRESS WITH COVER 195x85x18 cm - load 250 kg



Product Code	44803
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	2207518
CND	V080702
EAN/UPC	8023279448030

Description

MATTRESS WITH COVER 195 x 85 x 18 cm

LOAD UP TO 250 KG

Polyurethane foam mattress, density 35 kg/mc.
Non deformable, antibacterial, antifungal and autoclavable at 135°C.
Supplied with TEKMA polyurethane cover with zip fastener.