

KING PLUS EXAMINATION COUCH - blue



Product Code	44490
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No
RDM (NSIS)	2254329
CND	V0899
EAN/UPC	8023279444902

Description

KING PLUS EXAMINATION COUCH - blue

Strong and stylish structure in oval tube 50 x 25 mm painted with metallic grey epoxy powder. Top in washable, fire proof, leatherette material class 1 IM with high density padding. Head section adjustable by mechanical system, maximum 75°.

Mattress thickness: 7 cm with lateral white band.
Maximum load: 150 kg.
Size: 190 x 60 x h 80 cm
Delivered in kit form.

Multilingual instructions: GB, FR, IT, ES, PT, DE, NL, PL, RO, SE, GR, HR, SI, Arabic.

Made in Italy.