



Product Code	44900
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

EAN/UPC 8023279449006

## Description

S/S AISI 410 RUBBISH BIN  
Capacity: 8 l  
Bin size: Ø22.5 × h33 cm  
Bucket size: Ø21 × h27.5 cm  
Weight: 1.2 kg

Rubbish bins made of high-grade stainless-steel frame with mirror surface polishing, equipped with handle and foot pedal. Inner plastic bucket has strong pressure resistance, doesn't leak water and is easy to remove thanks to the handle. Thickness: 0.35 mm.

Inner bucket included.