

# ISO VISITOR CHAIR - fabric - red



|               |                   |
|---------------|-------------------|
| Product Code  | 45030             |
| Unit of sale  | 1 pc              |
| Minimum order | 1                 |
| Type          | No medical device |
| Class         |                   |

|         |               |
|---------|---------------|
| UK-REP  | Not required  |
| CH-REP  | Not required  |
| EAN/UPC | 8023279450309 |

## Description

Comfortable upholstered chair with 4 legged frame.

Can be used for patient and waiting rooms, with Cuneo, Cremona and Varese chairs.

Fire proof material.

Comfortable seat made of expanded polyurethane hit injected heliocell foams.

Maximum load: 110 kg.

Made in Italy.

## Technical Specifications

- Backrest height: 82 cm
- Seat height: 45 cm
- Width: 53 cm
- Depth: 43 cm
- Weight: 7 Kg