

ISO VISTOR CHAIR - fabric - blue



Product Code	45032
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279450323

Description

Comfortable upholstered chair with 4 legged frame.

Can be used for patient and waiting rooms, with Cuneo, Cremona and Varese chairs.

Fire proof material.

Comfortable seat made of expanded polyurethane hit injected heliocell foams.

Maximum load: 110 kg.

Made in Italy.

Technical Specifications

- Backrest height: 82 cm
- Seat height: 45 cm
- Width: 53 cm
- Depth: 43 cm
- Weight: 7 Kg