

CUNEO CHAIR without armrest - fabric- blue



Product Code	45060
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279450606

Description

Functional office upholstered chair, polypropylene swivel base, movable backrest with permanent contact.

Made in Italy.

- Fabric
- Blue

Technical Specifications

- Weight: 12 Kg
- Backrest height: 98/111 m
- Seat height: 42/52 cm
- Width: 56 cm
- Depth: 49 cm
- Maximum load: 110 kg