



Product Code	45065
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279450651

## Description

Functional office upholstered chair, with polypropylene swivel base, movable backrest with permanent contact.

- Leatherette
- Black

## Technical Specifications

- Weight: 12 Kg
- Backrest height: 98/111 m
- Seat height: 42/52 cm
- Width: 56 cm
- Depth: 49 cm
- Maximum load: 110 kg