

CUNEO CHAIR with armrest - leatherette - blue



Product Code 45076

Unit of sale 1 pc

Minimum order 1

Type No medical device

Class

UK-REP Not required

CH-REP Not required

EAN/UPC 8023279450767

Description

Functional office upholstered chair, with armrests and polypropylene swivel base, movable backrest with permanent contact.

Made in Italy.

- Leatherette
- Blue

Technical Specifications

- Weight: 12 Kg
- Backrest height: 98/111 m
- Seat height: 42/52 cm
- Width: 56 cm
- Depth: 49 cm
- Maximum load: 110 kg