

CREMONA CHAIR - fabric - grey



Product Code	45081
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279450811

Description

Superior quality chair with neat linens, upholstered seat and backrest for greater comfort, swivel base, armrests, permanent contact movement of the backrest.
High density foam. UNI EN 1335/1/2/3.

Maximum load: 110 kg

Made in Italy.

- Fabric
- Grey

Technical Specifications

- Weight: 12 Kg
- Backrest height: 90/107 cm
- Seat height: 45/52 cm
- Width: 58 cm
- Depth: 44 cm