



Product Code	45083
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279450835

## Description

Superior quality chair with neat linens, upholstered seat and backrest for greater comfort, swivel base, armrests, permanent contact movement of the backrest.  
High density foam. UNI EN 1335/1/2/3.

Maximum load: 110 kg

Made in Italy.

- Fabric
- Red

## Technical Specifications

- Weight: 12 Kg
- Backrest height: 90/107 cm
- Seat height: 45/52 cm
- Width: 58 cm
- Depth: 44 cm