

VARESE CHAIR with armrest - fabric - red



Product Code	45095
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279450958

Description

Quality chair with upholstered seat and backrest, permanent contact movement of backrest and swivel base.
Maximum load: 110 kg.

- Fabric
- Red

Made in Italy.

Technical Specifications

- Weight: 11.6 kg
- Backrest height: 94/112 cm
- Seat height: 46/59 cm
- Width: 60 cm
- Depth: 47 cm