

SIDNEY EXECUTIVE CHAIR



Product Code	45110
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

EAN/UPC	8023279451108
---------	---------------

Description

SIDNEY EXECUTIVE CHAIR

Stylish and comfort executive chairs with bracket to adjust height and back.
Interwoven mesh cloth, skin friendly, smooth, breathable and easy to clean with plastic armrests.
Chair base Ø 58 cm with chromed plated frame.
Five nylon 360° rotating castors.

Made in Italy.

Technical Specifications

Backrest size: 49x48 cm
Seat height: 45/55 cm
Width: 48 cm
Depth: 50 cm