

ELITE LOW-BACKED CHAIR - fabric - red



Product Code	45124
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279451245

Description

Executive low-backed armchair.
Steel and chrome plated base and arms.
Height adjustable by gas pump.

Made in Italy.

Technical Specifications

- Lining material: red fabric
- Weight: 14.5 kg
- Backrest height: 83/93 cm
- Seat height: 39/49 cm
- Width: 51 cm
- Depth: 45 cm
- Maximum load: 110 kg