



Product Code	45136
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279451368

## Description

Executive high-backed armchair.  
Steel and chrome plated base and arms.  
Height adjustable by gas pump.

Made in Italy.

## Technical Specifications

- Lining material: black leatherette
- Weight: 14.5 kg
- Backrest height: 101/111 cm
- Seat height: 39/49 cm
- Width: 51 cm
- Depth: 45 cm
- Maximum load: 110 kg