

SURGILANCE SAFETY NEEDLE G21 - automatic lancets



Product Code	24517
Unit of sale	box of 100 pcs
Minimum order	1
Type	Medical device
Class	II A

UK-REP	Yes
CH-REP	No
RDM (NSIS)	16427
CND	V9009
EAN/UPC	8023279245172
UMDNS	16380
PARAF	977810682
GMDN	61578

Description

SURGILANCE™ SAFETY AUTOMATIC LANCETS - NEEDLES

Colour: Pink
Depth of penetration: 2.8 mm
Entire needle Gauge: 21G
Needle tip Gauge: 21G
Blood Flow: 40-60 µl
Finger Stick: yes
Glucose: yes
Hematocrit: yes

Multi purpose capillary blood sampling devices for all testing requirements.
The one-step safety lancet is easier to use, no arming is required, and safer, once lancet is used it is rendered inoperable.
High-speed delivery and penetration method minimizes patient pain and operator error.
Latex-free, hypo-allergenic, disposable, sterile.

Multilingual box: GB, FR, IT, ES, DE, PT, NL, SE, DK, RU, GR