



| | |
|---------------|-------------------|
| Product Code | 45170 |
| Unit of sale | 1 pc |
| Minimum order | 1 |
| Type | No medical device |
| Class | |

| | |
|---------|---------------|
| UK-REP | Not required |
| CH-REP | Not required |
| EAN/UPC | 8023279451702 |
| PARAF | 977816913 |

Description

ARIANNA LIFT ARMCHAIRS WITH 1 MOTOR

1 motor lift armchair: the chair's backrest and footrest can lift or recline together.

- Color: blue
- 1 motor
- Fireproof: NO
- Footrest: standard

Lift armchairs that feature a powered lifting mechanism that pushes the chair up from its base to assist the user to a standing position.

Power supply: 100-240 V - 50/60 Hz.

Multilingual manual: GB, FR, IT, ES.

Technical Specifications

- Seat
Width: 49 cm
Depth: 49.5 cm
Height: 49 cm
- Backrest height: 70 cm
- Chair
Width: 75 cm
Length (reclined): 170.5 cm
- Load: 150 Kg (331 lbs)