



Product Code	45177
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279451771
PARAF	977816952

## Description

### ARIANNA LIFT ARMCHAIRS WITH 2 MOTORS

2 motors lift armchair with standard or extra lifting footrest: it can reach the bed position providing a relaxing experience. Its backrest and footrest can operate separately to any position. Extra lifting footrest allows 8 cm lifting adjustment higher than flat position.

- Color: green
- 2 motors
- Fireproof: NO
- Footrest: standard

Lift armchairs that feature a powered lifting mechanism that pushes the chair up from its base to assist the user to a standing position.

Power supply: 100-240 V - 50/60 Hz.

Multilingual manual: GB, FR, IT, ES.

## Technical Specifications

- Seat  
Width: 49 cm  
Depth: 49.5 cm  
Height: 49 cm
- Backrest height: 70 cm
- Chair  
Width: 75 cm  
Length (reclined): 170.5 cm
- Load: 150 Kg (331 lbs)