

GINEVRA LIFT ARMCHAIR 2 motors - blue



Product Code	45180
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279451801
PARAF	977816976

Description

GINEVRA LIFT ARMCHAIRS WITH 2 MOTORS
2 motors lift armchair with standard or extra lifting footrest: it can reach the bed position providing a relaxing experience. Its backrest and footrest can operate separately to any position. Extra lifting footrest allows 8 cm lifting adjustment higher than flat position.

- Color: blue
- 2 motors
- Fireproof: NO
- Footrest: standard

Lift armchairs that feature a powered lifting mechanism that pushes the chair up from its base to assist the user to a standing position.

Power supply: 100-240 V - 50/60 Hz.

Multilingual manual: GB, FR, IT, ES.

Technical Specifications

- Seat
Width: 46 cm
Depth: 51 cm
Height: 40 cm
- Backrest height: 66 cm
- Chair
Width: 74 cm
Lenght (reclined): 169 cm
- Load: 150 Kg (331 lbs)